

Table of contents

\\USER

Yamin, Stephanie

SYA-COMPAS

CDIPv3.8

12Ch	multiplane	localizer
3D_T1		
PD_T2		
2D_FLAIR_FS		
T2-star		
ep2d_diff_30dir_AP		
ep2d_diff_PA		
ep2d_bold_resting_state		
FieldMap		

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\12Ch multiplane localizer

TA: 0:11 PM: REF Voxel size: 1.3×1.3×7.5 mmPAT: 2 Rel. SNR: 1.00 : tfl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	9
Dist. factor	42 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	9
Dist. factor	42 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice group	3
Slices	9
Dist. factor	42 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	7.5 mm
TR	262.99 ms
TE	1.41 ms
Averages	1
Concatenations	27
Filter	Prescan Normalize
Coil elements	HEP;NEP

Contrast - Common

TR	262.99 ms
TE	1.41 ms
TD	0 ms
Magn. preparation	None
Flip angle	50 deg
Fat suppr.	None
Wrap-up Magn.	Restore

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	280 mm
FoV phase	100.0 %
Slice thickness	7.5 mm
Base resolution	208
Phase resolution	66 %
Phase partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	46
Matrix Coil Mode	Auto (Triple)
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slice group	1
Slices	9
Dist. factor	42 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	9
Dist. factor	42 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice group	3
Slices	9
Dist. factor	42 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	7.5 mm
TR	262.99 ms
Multi-slice mode	Sequential
Series	Descending
Concatenations	27

Geometry - AutoAlign

Slice group	1
Slice group	2
Slice group	3

Geometry - AutoAlign

AutoAlign	---
Position	L0.0 A30.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Initial Position	L0.0 A30.0 H0.0
L	0.0 mm
A	30.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	Restore
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Coil Mode	Auto (Triple)
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	On
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	262.99 ms
Concatenations	27
Segments	69

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	280 mm
FoV phase	100.0 %
Phase resolution	66 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	0

Physio - PACE

Resp. control	Off
Concatenations	27

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	1
TE	1.41 ms
TR	262.99 ms
Save Original Images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save Original Images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Echo spacing	3.2 ms
Sequence type	Trufi
Bandwidth	1093 Hz/Px

Sequence - Part 2

Define	Shots
Shots per slice	1
Segments	69
Trufi delta freq.	0 Hz

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
Cine	Off

Sequence - Assistant

Mode	Min flip angle
Min flip angle	45 deg
Allowed delay	5 s

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\3D_T1

TA: 5:21 PM: FIX Voxel size: 1.0×1.0×1.0 mmPAT: 2 Rel. SNR: 1.00 : tfl

Properties

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	192
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2300.0 ms
TE	2.98 ms
Averages	1
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP;NEP

Contrast - Common

TR	2300.0 ms
TE	2.98 ms
Magn. preparation	Non-sel. IR
TI	900 ms
Flip angle	9 deg
Fat suppr.	None
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off

Resolution - Common

Interpolation	Off
---------------	-----

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	192
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2300.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
AutoAlign	Head > Brain
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Coil Mode	Auto (Triple)
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
A >> P	256 mm
F >> H	256 mm
R >> L	192 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	2300.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	Non-sel. IR
TI	900 ms
Fat suppr.	None
Dark blood	Off
FoV read	256 mm
FoV phase	100.0 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save Original Images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save Original Images	On
MapIt	None
Flip angle	9 deg
Measurements	1
TR	2300.0 ms
TE	2.98 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	7.1 ms
Bandwidth	240 Hz/Px

Sequence - Part 2

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	192

Sequence - Assistant

Mode	Off
------	-----

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\PD_T2

TA: 3:11 PM: REF Voxel size: 0.9×0.9×3.0 mmPAT: 2 Rel. SNR: 1.00 : tse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	3000.0 ms
TE 1	10 ms
TE 2	93 ms
Averages	1
Concatenations	3
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP;NEP

Contrast - Common

TR	3000.0 ms
TE 1	10 ms
TE 2	93 ms
TD	0.0 ms
MTC	Off
Magn. preparation	None
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Water suppr.	None
Restore magn.	Off

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
Base resolution	256
Phase resolution	100 %

Resolution - Common

Phase partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	On

Geometry - Common

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	R >> L
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	3000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	3

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	R >> L
Initial Position	L0.7 A14.6 F9.0
L	0.7 mm
A	14.6 mm
F	9.0 mm
Initial Rotation	90.00 deg
Initial Orientation	T > C
T > C	-9.0
> S	0.1

Geometry - Saturation

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Water suppr.	None
Restore magn.	Off
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

Geometry - Tim CT

Tim CT mode	Off
Slices	48
Slice thickness	3.0 mm
Dist. factor	0 %
FoV read	240 mm
FoV phase	100.0 %

System - Miscellaneous

Positioning mode	REF
Table position	F
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Coil Mode	Auto (Triple)
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Rotation	90.00 deg
R >> L	240 mm
A >> P	240 mm
F >> H	144 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	3000.0 ms
Concatenations	3

Physio - Cardiac

Magn. preparation	None
Fat suppr.	Fat sat.

Physio - Cardiac

Dark blood	Off
FoV read	240 mm
FoV phase	100.0 %
Phase resolution	100 %
Trajectory	Cartesian

Physio - PACE

Resp. control	Off
Concatenations	3

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save Original Images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	Off
Contrasts	2
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	10.3 ms
Bandwidth	181 Hz/Px

Sequence - Part 2

Define	Turbo factor
Echo trains per slice	20
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast
Hyperecho	Off
WARP	Off
Red. EC sensitivity	Off
Turbo factor	7

Sequence - Assistant

Mode	Min flip angle
Min flip angle	130 deg
Allowed delay	60 s

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\2D_FLAIR_FS

TA: 2:44 PM: FIX Voxel size: 0.9×0.9×3.0 mmPAT: 2 Rel. SNR: 1.00 : tir

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	9000.0 ms
TE	120.0 ms
Averages	1
Concatenations	2
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP;NEP

Contrast - Common

TR	9000.0 ms
TE	120.0 ms
TD	0.0 ms
MTC	Off
Magn. preparation	Slice-sel. IR
TI	2500 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Water suppr.	None
Restore magn.	Off
Freeze suppressed tissue	Off

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
Base resolution	256
Phase resolution	100 %

Resolution - Common

Phase partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	47
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	On

Geometry - Common

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	R >> L
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	9000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	2

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	R >> L
Initial Position	L0.7 A14.6 F9.0
L	0.7 mm
A	14.6 mm
F	9.0 mm
Initial Rotation	90.00 deg
Initial Orientation	T > C
T > C	-9.0
> S	0.1

Geometry - Saturation

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Water suppr.	None
Restore magn.	Off
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

Geometry - Tim CT

Tim CT mode	Off
Slices	48
Slice thickness	3.0 mm
Dist. factor	0 %
FoV read	240 mm
FoV phase	100.0 %

System - Miscellaneous

Positioning mode	FIX
Table position	F
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Coil Mode	Auto (Triple)
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Rotation	90.00 deg
R >> L	240 mm
A >> P	240 mm
F >> H	144 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	9000.0 ms
Concatenations	2

Physio - Cardiac

Magn. preparation	Slice-sel. IR
TI	2500 ms

Physio - Cardiac

Fat suppr.	Fat sat.
Dark blood	Off
FoV read	240 mm
FoV phase	100.0 %
Phase resolution	100 %
Trajectory	Cartesian

Physio - PACE

Resp. control	Off
Concatenations	2

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save Original Images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	9.26 ms
Bandwidth	222 Hz/Px

Sequence - Part 2

Define	Turbo factor
Echo trains per slice	8
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast
Hypercho	Off
WARP	Off
Red. EC sensitivity	Off
Turbo factor	19

Sequence - Assistant

Mode	Off
Allowed delay	60 s

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\T2-star

TA: 3:04 PM: FIX Voxel size: 0.9×0.9×3.0 mmPAT: 2 Rel. SNR: 1.00 : fl

Properties

Prio recon	On
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Routine

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	650.0 ms
TE	20.00 ms
Averages	1
Concatenations	2
Filter	None
Coil elements	HEA;HEP;NEP

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	650.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	2

Contrast - Common

TR	650.0 ms
TE	20.00 ms
TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
Initial Position	L0.7 A14.6 F9.0
L	0.7 mm
A	14.6 mm
F	9.0 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-9.0
> S	0.1

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

Geometry - Tim CT

Tim CT mode	Off
Slices	48
Slice thickness	3.0 mm
Dist. factor	0 %
FoV read	240 mm
FoV phase	100.0 %
Segments	1

System - Miscellaneous

Positioning mode	FIX
Table position	F
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Coil Mode	Auto (Triple)
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	650.0 ms
Concatenations	2
Segments	1

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	240 mm
FoV phase	100.0 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	2

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save Original Images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save Original Images	On

Inline - Soft Tissue

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save Original Images	On
MapIt	None
Flip angle	20 deg
Measurements	1
Contrasts	1
TR	650.0 ms
TE	20.00 ms

Sequence - Part 1

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Bandwidth	200 Hz/Px

Sequence - Part 2

Segments	1
Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\ep2d_diff_30dir_AP

TA: 5:27 PM: FIX Voxel size: 2.0×2.0×2.0 mmPAT: 2 Rel. SNR: 1.00 : epse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	70
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	9200 ms
TE	93.0 ms
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP;NEP

Contrast - Common

TR	9200 ms
TE	93.0 ms
MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Contrast - Dynamic

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36

Resolution - iPAT

Matrix Coil Mode	Auto (Triple)
Reference scan mode	GRE/separate

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	On
Dynamic Field Corr.	Off

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	70
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	9200 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
Initial Position	L0.7 A14.6 F9.0
L	0.7 mm
A	14.6 mm
F	9.0 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-9.0
> S	0.1

Geometry - Saturation

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm

System - Miscellaneous

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Matrix Coil Mode	Auto (Triple)
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Rotation	0.00 deg
A >> P	256 mm
R >> L	256 mm
F >> H	140 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	9200 ms
Concatenations	1

Physio - PACE

Resp. control	Off
Concatenations	1

Diff - Neuro

Diffusion mode	MDDW
Diff. directions	30
Diffusion Scheme	Monopolar
Diff. weightings	2
b-value 1	0 s/mm ²
b-value 2	1000 s/mm ²
b-value 1	3
b-value 2	1
Diff. weighted images	On
Trace weighted images	On
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Noise level	40

Diff - Body

Diffusion mode	MDDW
----------------	------

Diff - Body

Diff. directions	30
Diffusion Scheme	Monopolar
Diff. weightings	2
b-value 1	0 s/mm ²
b-value 2	1000 s/mm ²
b-value 1	3
b-value 2	1
Diff. weighted images	On
Trace weighted images	On
ADC maps	Off
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	Off
b-Value >=	0 s/mm ²
Noise level	40

Diff - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Optimization	None
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.65 ms
Bandwidth	1776 Hz/Px

Sequence - Part 2

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\ep2d_diff_PA

TA: 1:09 PM: FIX Voxel size: 2.0×2.0×2.0 mmPAT: 2 Rel. SNR: 1.00 : epse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	70
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0
Phase enc. dir.	P >> A
AutoAlign	---
Phase oversampling	0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	9200 ms
TE	93.0 ms
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP;NEP

Contrast - Common

TR	9200 ms
TE	93.0 ms
MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Contrast - Dynamic

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36

Resolution - iPAT

Matrix Coil Mode	Auto (Triple)
Reference scan mode	GRE/separate

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	On
Dynamic Field Corr.	Off

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	70
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0
Phase enc. dir.	P >> A
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	9200 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0
Phase enc. dir.	P >> A
Initial Position	L0.7 A14.6 F9.0
L	0.7 mm
A	14.6 mm
F	9.0 mm
Initial Rotation	-180.00 deg
Initial Orientation	T > C
T > C	-9.0
> S	0.0

Geometry - Saturation

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm

System - Miscellaneous

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Matrix Coil Mode	Auto (Triple)
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0
Rotation	180.00 deg
A >> P	256 mm
R >> L	256 mm
F >> H	140 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	9200 ms
Concatenations	1

Physio - PACE

Resp. control	Off
Concatenations	1

Diff - Neuro

Diffusion mode	MDDW
Diff. directions	30
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	0 s/mm ²
b-value	5
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	40

Diff - Body

Diffusion mode	MDDW
Diff. directions	30
Diffusion Scheme	Monopolar

Diff - Body

Diff. weightings	1
b-value	0 s/mm ²
b-value	5
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	Off
b-Value >=	0 s/mm ²
Noise level	40

Diff - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Optimization	None
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.65 ms
Bandwidth	1776 Hz/Px

Sequence - Part 2

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\ep2d_bold_resting_state

TA: 9:03 PM: FIX Voxel size: 3.5×3.5×3.5 mmPAT: 2 Rel. SNR: 1.00 : epfid

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	40
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
TR	2140 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP;NEP

Contrast - Common

TR	2140 ms
TE	30.0 ms
MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	250
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2

Resolution - iPAT

Ref. lines PE	24
Matrix Coil Mode	Auto (Triple)
Reference scan mode	EPI/separate

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	On

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry - Common

Slice group	1
Slices	40
Dist. factor	0 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
TR	2140 ms
Multi-slice mode	Interleaved
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
Initial Position	L0.7 A14.6 F9.0
L	0.7 mm
A	14.6 mm
F	9.0 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-9.0
> S	0.1

Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

System - Miscellaneous

Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Coil Mode	Auto (Triple)
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

Sequence - Part 2

EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

System - Adjustments

B0 Shim mode	Advanced
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Rotation	0.00 deg
A >> P	224 mm
R >> L	224 mm
F >> H	140 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	2140 ms
Concatenations	1

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial filter	Off
Measurements	250
Delay in TR	0 ms
Multiple series	Off

Sequence - Part 1

Introduction	On
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.5 ms
Bandwidth	2442 Hz/Px

\\USER\Yamin, Stephanie\SYA-COMPAS\CDIPv3.8\FieldMap

TA: 1:13 PM: FIX Voxel size: 3.0×3.0×3.5 mmRel. SNR: 1.00 : fm_r

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	40
Dist. factor	10 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
TR	476.0 ms
TE 1	4.84 ms
TE 2	7.3 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HEA;HEP;NEP

Contrast - Common

TR	476.0 ms
TE 1	4.84 ms
TE 2	7.3 ms
MTC	Off
Flip angle	50 deg
Fat suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
Base resolution	74
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

Matrix Coil Mode	Auto (Triple)
------------------	---------------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	40
Dist. factor	10 %
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
TR	476.0 ms
Multi-slice mode	Interleaved
Series	Descending
Concatenations	1

Geometry - AutoAlign

Slice group	1
AutoAlign	---
Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Phase enc. dir.	A >> P
Initial Position	L0.7 A14.6 F9.0
L	0.7 mm
A	14.6 mm
F	9.0 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-9.0
> S	0.1

Geometry - Saturation

Fat suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off

System - Miscellaneous

Matrix Coil Mode	Auto (Triple)
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	L0.7 A14.6 F9.0 mm
Orientation	T > C-9.0 > S0.1
Rotation	0.00 deg
A >> P	224 mm
R >> L	224 mm
F >> H	154 mm
Reset	Off

System - Tx/Rx

Frequency 1H	123.259052 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Introduction	On
Dimension	2D
Asymmetric echo	Off
Contrasts	2
Flow comp.	Yes
Multi-slice mode	Interleaved
Bandwidth	295 Hz/Px

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

Sequence - Assistant

Mode	Off
------	-----